


# Mini VACUUM PUMP



### MINI VACUUM PUMP LMVP-A10

Model of Labtron Mini Vacuum Pump LMVP-A Series uses big preface torque configuration and is particularly designed for low temperature working environment. It portable and requires less storage place.

#### Features

- ▶ Compact, reliable and powerful
- ▶ Especially designed for low temperature working environment
- ▶ Space-saving and ideal for systems integration
- ▶ High pumping speeds with small dimensions
- ▶ Eco-friendly exhaust gas free of oil mists
- ▶ Quiet, low-vibration operation
- ▶ Easy to maintain

#### Application

It can be used for vacuum packaging, vacuum freezing, gas analysis and other vacuum evacuations.



LMVP-A10

## Specification

Model	LMVP-A10
Pumping Speed	28.31 to 22.65L/min
Ultimate Pressure	5Pa
Weight	6.5kg
Rotating Speed	1720/1440rpm
Oil Capacity	220ml
Dimension	243x114x207mm
Power	90/120W
Power Supply	110V, 60Hz/220V, 50Hz



---

**LABTRON EQUIPMENT LTD**

Basing View, Basingstoke, Hampshire RG21 4RG  
UK

 01256 610 110 |  [info@Labtron.com](mailto:info@Labtron.com) |  [www.Labtron.com](http://www.Labtron.com)